







# One Step Pregnancy Midstream Test

Start testing from 7 days past ovulation!

Fertility2Family One Step Pregnancy Midstream Tests are a fast, reliable and easy-to-use method of detecting pregnancy. Fertility2Family Pregnancy Midstream Test detects the hCG (human Chorionic Gonadotropin) hormone. A hormone which is produced by the placenta after the fertilized egg adheres to the wall of the uterus.

## How to Perform a Pregnancy Test?

1. Remove the test stick from the foil wrapper.
2. Remove the cap to expose the absorbent tip, avoiding contact with the absorbent end.
3. Hold the stick by the thumb grip with the exposed absorbent tip pointing downward.  
Urinate onto the absorbent tip only until it is thoroughly wet.
4. Recap the device and lay the test stick down on a flat surface with the window facing up.
5. Results may be observed after 1 minute. Do not interpret test results after 5 minutes.
6. Discard the test device after single use.



## Interpreting the Test Results

### Negative:

If only one coloured line appears in the control area, the test is negative, and you can assume that you are not pregnant.

### Positive:

If two coloured lines appear both in the control area and test area, the test is positive, and you can assume you are pregnant.

### Inconclusive:

If there is no distinct pink-purple colored band visible both in the test area and control area, the test is inconclusive. It is recommended that you should retest, using a new test stick. If the problem persists, contact Fertility2Family.

## FAQs

### Q: When can I start testing?

A: There is no need to wait until you have missed your period to begin testing. You can test as early as 7 days past ovulation, or 5-7days before your expected period. When testing early, try to use first morning urine.

### Q: What should I do if I get a positive result?

A: You should make an appointment to see your doctor to confirm the pregnancy and begin prenatal care.

### Q: What should I do if I get a negative result?

A: A negative result indicated that the hCG hormone was not detected in your urine, and you are unlikely pregnant. If your period is late, we suggest you re-test, and if this test is also negative we recommend you contact your doctor.

### Q: The Test Line is not as dark as the Control Line. Am I pregnant?

A: Yes, you are pregnant even if a very light test line is detected.

### Q: What can I do to increase my chances of falling pregnant?

A: Primarily you need to determine when you are ovulating, as this is your most fertile window. By using Fertility2Family **Ovulation Tests**, or our **Basal Body Temperature(BBT) thermometer** you can determine your most fertile window, and dramatically increase your chances of falling pregnant. You may also find our **Sperm Friendly Lubricant** helpful. This lubricant has been scientifically formulated to match fertile cervical liquid and sperm, enhance sexual activity and increase the likelihood of conception.

You can purchase these products at [www.fertility2family.com.au](http://www.fertility2family.com.au)

